

## HANDBOOK

# Bio-T kit® *Dictyocaulus viviparus*

Cat. N° BIOTK096 - 50 reactions

**Detection of *Dictyocaulus viviparus* (*D. viviparus*)  
by real-time PCR (qPCR)  
with Exogenous internal positive control (IPC)**

## BOVINE

### Sample types

- Faeces
- Individual analysis or by pool up to 5

### Recommended nucleic acid (NA) extractions

- Silica membrane columns extraction (eg : BioSellal – BioExtract® Column Cat. N° BEC050 or BEC250)
- Magnetic beads extraction (eg : BioSellal – BioExtract® SuperBall® Cat. N° BES384)

*Veterinary use only*



## DOCUMENTS MANAGEMENT

The Bio-T kit® *Dictyocaulus viviparus* has two technical handbooks:

- The extraction handbook shared between all the Bio-T kit® of the STOOL line, displaying BioSella's recommended extraction protocols.
- The Bio-T kit® *Dictyocaulus viviparus* qPCR handbook, presenting the instruction information to perform the qPCR.

The last versions in use for each handbook are indicated on the certificate of analysis (CA) provided with the Bio-T kit® *Dictyocaulus viviparus*.

Besides these two handbooks, a summary report of the validation file is available on request, contact BioSella (contact@biosella.com).

## MODIFICATIONS MANAGEMENT

BioSella indicates modifications done to this document by highlighting them using the rules presented in the Table below:

MODIFICATIONS MANAGEMENT			
Type of modification Highlighting color	Minor modifications	Type 1 Major modifications	Type 2 Major modifications
Impact on revision / version	Change of revision date No change of version	Change of revision date + change of version	Change of revision date + change of version
Examples of modifications	Corrections: typographical, grammatical or turns of phrase	EPC reference modification	Modification of Master Mix composition
	Addition of new sample type for extraction	Exogenous IPC reference modification	Modification of validated extraction protocol
	Addition of information giving more details or alternative protocol		
	Addition/Suppression of optional information		

## PRESENTATION

### Recommendations for sampling, shipping and storage of samples

Real-time PCR is a powerful technique allowing the detection of few amounts of pathogen genome. Genome can be rapidly degraded depending on the pathogen nature (bacteria / parasites, enveloped viruses...), the genome nature (DNA / RNA) and the sample type (presence of DNase / RNase). Thus, BioSellal recommends the following instructions to guarantee an optimal diagnosis.

#### Sampling

To prevent cross-contamination between samples leading to false positive results, it is mandatory to use disposable materials for single use and to avoid direct contact between specimens.

#### Shipping

It is mandatory to ship within 7 days after collection. The shipment must be done, if possible, under cover of the positive cold.

#### Storage after reception

Recommended storage of samples at  $5^{\circ}\text{C} \pm 3$  for a maximum of 7 days and  $\leq -16^{\circ}\text{C}$  beyond.

## STOOL Line

This kit belongs to the STOOL line which gather a set of kits sharing common extraction and qPCR protocols. It is compatible with BioSellal's other kits except with the ones belonging to the PIG and AVIAN lines. (information available on [www.biosellal.com](http://www.biosellal.com)).

## Description of the Bio-T kit® *Dictyocaulus viviparus*


The **Bio-T kit® *Dictyocaulus viviparus*** (Cat. N° BIOTK096) contains a ready to use **PCR Master Mix** allowing the detection **in the same reaction well of**:

- ***Dictyocaulus viviparus* (D. viviparus)** with a 6-FAM labelling
- An **Exogenous internal positive control IPC DNA**, with a Cy5 labelling, to add at the extraction step to assess nucleic acids extraction quality and absence of PCR inhibitors.

This kit, based on qualitative detection (detected or non detected) from Faeces samples (Individual analysis or by pool up to 3), was developed and validated according to the **French regulatory standard NF U47-600-2 edited by AFNOR** for the PCR part.

**Extraction protocols recommended by BioSella are described in in the extraction handbook of the STOOL line.**

## Description of the whole process



STOOL line Extraction handbook		qPCR handbook of the Bio-T kit® <i>Dictyocaulus viviparus</i>		
Faeces *	BioExtract® Column BioExtract® SuperBall®	Ready-to-use Master Mix MMDvivip-A	Samples NC/NCS Process positive control EPC (EPCDvivip-A)	Dyes: FAM/Cy5 Passive reference: ROX Program: Classical program ± RT Standard or Fast ramping

\* pretreatment mandatory

## Kit contents and storage

**Table 1. Description of the kit contents**

Description	Reference	Volume/tube	Presentation	Storage
<b>Master Mix (MM)</b> Ready to use	MMDvivip-A	750 µl	<b>White</b> cap tube Bag A	≤-16°C Protected from light, « MIX » Area
<b>Internal Positive Control (IPC)</b> Exogenous Internal Positive Control	IPC-A	250 µl	<b>Pink</b> cap tube Bag B	≤-16°C « Extraction » Area
<b>External Positive Control (EPC)</b> Positive PCR control of <i>D. viviparus</i>	EPDVDvivip-A	110 µl	<b>Orange</b> cap tube Bag C	≤-16°C « Addition of Nucleic acids » Area
<b>Water</b> RNase/DNase free	Aqua-A	1 ml	<b>Blue</b> cap tube Bag C	≤-16°C « Addition of Nucleic acids » Area

## List of consumables and reagents not included in kit

Table 2. Consumables and reagents not included in kit			
Consumables/ Reagents	Description	Provider	Cat. N°
ATL Buffer	Lysis Buffer	BioSellal	ATL19076
BioExtract® Column	DNA/RNA column extraction kit (50)	BioSellal	BEC050
BioExtract® Column	DNA/RNA column extraction kit (250)	BioSellal	BEC250
BioExtract® SuperBall®	DNA/RNA Magnetic beads extraction kit (4 x 96)	BioSellal	BES384
BioPrep Parasito - A	Sample prep extraction kit conical bottom tube (50) Conical bottom tube 50 ml – BioPrepParasito - 1 Tube 2 ml – BioPrepParasito - 2 45 ml ATL Buffer	BioSellal	BPPARASIT01A
BioPrep Parasito - B	Sample prep extraction kit flat bottom tube (50) Flat bottom tube 50ml – BioPrepParasito - 1 Tube 2ml – BioPrepParasito - 2 45 ml ATL Buffer	BioSellal	BPPARASIT01B

For consumables related to the thermal cycler, refer to the user manual of the device.

## List of reagents for performances validation

To confirm performances of your thermocycler(s), *D. viviparus* DNA (quantified in copies/qPCR) provided with the qPCR kit (orange cap tube) could be used. Please, contact BioSellal for more information (tech@biosellal.com).

## Main critical points

- Wear appropriate personal protective equipment (lab coat, disposable gloves frequently changed).
- Work in dedicated and separate areas to avoid contamination: "Extraction" (unextracted samples storage, extraction equipment area), "Mix" (ready to use MM storage, qRT-PCR plates preparation), "Nucleic acids (NA) Addition" (Nucleic Acid storage and addition of extracted NA and controls in the qRT-PCR plate), "PCR" (final area containing the thermocycler(s)).
- Use dedicated equipment for each working area (gloves, lab coat, pipettes, vortex, ...).
- Use filter tips.
- Before use, thaw all components at room temperature.
- Vortex and spin briefly (mini-centrifuge) all reagents before use.
- Avoid the repetition of freezing-thawing cycles for samples, lysates, extracted NA.
- **Pathogen's genome detected by the STOOL line's kits can be DNA or RNA. Working with RNA is more demanding than working with DNA** (RNA instability and omnipresence of the RNases). For these reasons, special precautions must be taken:
  - o Always wear gloves, change them frequently, especially after contact with skin or work surfaces.
  - o Treat all surfaces and equipment with RNases inactivation agents (available commercially).
  - o When wearing gloves and after material decontamination, minimize the contact with surfaces and equipment in order to avoid the reintroduction of RNases.
  - o Use "RNase free" consumable.
  - o It is recommended to store the RNA at  $\leq 5 \pm 3^{\circ}\text{C}$  during the manipulation and then freeze it as soon as possible, preferably at  $\leq -65^{\circ}\text{C}$  or by default at  $\leq -16^{\circ}\text{C}$ .
  - o Open and close tubes one by one in order to limit the opening times and avoid any contact with RNases present in the environment (skin, dust, working surfaces...).

# DETECTION OF *D. VIVIPARUS* BY qPCR WITH BIOTK096

## Global Procedure

### 1) Establish qPCR plate setup defining each sample position and including the following controls:

- **Negative Control Sample (NCS):** water (or PBS) replaces the sample from the first step of sample preparation.


This control is mandatory for each extraction series.

- **Negative Amplification Control (NC):** 5 µl of water RNase/DNase free (Aqua-A, **blue** cap tube) replaces sample Nucleic Acid extract on qPCR plate.

This control is recommended when using the kit for the first time or to verify the absence of Master Mix contamination.

- **External Positive Control of *D. viviparus* (EPC):** Synthetic DNA (EPCDvivip-A, **orange** cap tube), containing specific target of *D. viviparus*

This control is mandatory.

 **CAUTION:** *EPC tube handling represents a nucleic acid contamination hazard, it is thus recommended to open and handle it in a restricted area, away from other PCR components and to take precautions to avoid cross-contamination with nucleic acids extracts during deposit on the qPCR plate.*

- If available, a **Process Positive Control (MRI)**, a weak positive sample of Faeces is extracted in parallel with tested samples. After qPCR, MRI Ct values will be monitored on a Shewhart control card. Obtaining conform Ct values validates the whole process. In this case, the use of the EPC, provided with the kit, is not mandatory.

## 2) qPCR plate preparation

### In the “MIX” dedicated area

1. After thawing, vortex and rapid centrifugation, **transfer 15 µl Master Mix MMDvivip-A (White cap)** in each well of interest (samples and controls).

### In the “Nucleic Acid addition” dedicated area

2. **Add 5 µl of extracted nucleic acids (or NCS, water, MRSI or EPC: EPCDvivip-A orange cap tube)** in each well of interest. Make sure to pipet out in the bottom of the well, in the Master Mix, and to avoid the formation of bubbles.

*Note:* if the exogenous IPC was not added during sample extraction, it can be added directly in the qPCR plate:

- Add 1 µl of IPC (**pink** cap) with the extracted nucleic acids
  - Or add directly the IPC (1 µl per reaction) in an aliquot of Master Mix before the deposits of 16 µl of this mix into each well of interest. Then add 5 µl of extracted nucleic acids.
- The reaction volume will be increased to 21 µl, without impacting the performances of the qPCR.

3. **Seal the plate with an optically clear sealer or close the strip caps.**

### In the “PCR” amplification dedicated area

4. **Define the thermocycler parameters** (see Table 3, Table 4, Table 5, Table 6)
5. It is recommended to **spin the plate down prior to place it in the thermocycler**, to prevent drops in the well pit walls.
6. Start the qPCR program. Approximate run time: 60 min.

### 3) Thermocycler settings

This kit was developed and validated on AriaMx™ (Agilent Technologies, Fast ramping by default) and confirmed on ABI PRISM® 7500 Fast (Applied Biosystems) in standard ramping and fast ramping. It is compatible with all thermal cyclers with at least 6-FAM and Cy5 channels. For more information, contact our technical support.

Table 3. Thermal cycler configuration		
	ABI PRISM® 7500 Fast	AriaMx™
<b>Mode</b>	Quantitation – Standard curve	Quantitative PCR, Fluorescence Probe
<b>Ramping</b>	Standard or Fast Ramping	Fast Ramping by default
<b>Passive Reference</b>	ROX	ROX

Table 4. Thermal cycler Settings			
Target	Detectors		Final Volume / well
	Reporter	Quencher	
<i>D. viviparus</i>	FAM	NFQ-MGB ou None*	20 µl
Exogenous IPC	Cy5	NFQ-MGB ou None*	= 15 µl Master Mix + 5 µl extracted nucleics acids or controls <sup>†</sup>
To assign to samples and controls <sup>†</sup>			

\* Depends on the thermocycler model. Do not hesitate to contact the BioSella Technical Support (tech@biosellal.com)

† Controls are NC (water), NCS (extracted water), EPC or extracted MRI.

Table 5. CLASSICAL Amplification program settings without RT		
Standard or Fast Ramping		
Cycles	Time	Temperature
1 cycle	5 min	95°C
40 cycles	15 sec	95°C
	30 sec* + data acquisition	60°C

\* Set 31s for some thermocyclers such as ABI PRISM® 7500.

† optional step, in case of simultaneous detection of RNA genomes. Achieving a reverse-transcription (RT) step prior to PCR for the amplification of RNA genomes has no impact on the performances of the Bio-T kit® *Dictyocaulus viviparus* (see the summary of the validation file).

NB: This amplification program is compatible with all Bio-T kit® except for ones belonging to the PIG and AVIAN LINES.

For thermal cycler such as LightCycler®480 and LightCycler®96 (Roche Life Science), it is recommended to use the following program:

Table 6. PIG/AVIAN Amplification program settings without RT <sup>†</sup>		
Standard Ramping		
Cycles	Time	Temperature
1 cycle	5 min	95°C
40 cycles	10 sec	95°C
	45 sec + data acquisition	60°C

<sup>†</sup> A reverse-transcription (RT) step prior to PCR can be performed in the context of a simultaneous analysis with pathogens whose genome is RNA.

NB: This amplification program is compatible with all Bio-T kit® of the PIG and AVIAN LINES.

## RESULTS INTERPRETATION

To analyze and interpret the signals obtained by qPCR, the Threshold must be set up.

The threshold must be assigned carefully in order to obtain the most reproducible result between different manipulations according to the requirements defined in Annex C of the French Standard **NF U47-600 (part 1)**. A consistent set of positives controls, usually an In-house Reference Material (MRI) or the EPC, is used to set the threshold value above the baseline and in the exponential amplification phase of the plot.

The Threshold Cycle, named « Ct » or « Cq » (depending on thermocyclers), corresponds to the intersection between the amplification curves and the threshold line. It allows the relative measurement of the concentration of the target in the PCR reaction when a calibrated extract is analyzed in the same series.

The qPCR series is validated if the controls (EPC, MRI, NCS and NC) present valid results, then the result of each sample can be interpreted.

## Main Scenarios

### Controls Reading

Table 7. PCR Controls results interpretation			
	Targets		Interpretation
	<i>D. viviparus</i> (FAM)	Exogenous IPC (Cy5)	
<b>NCS</b> Negative Control Sample  <b>MANDATORY</b>	Neg	Pos	Valid
	Pos	Pos	Contamination with a positive/negative sample during extraction step or during qPCR plate preparation.
	Neg	Neg	Omission of exogenous IPC addition? Defective extraction?
<b>NC</b> Negative PCR Control  <b>OPTIONAL</b>	Neg	Neg	Valid
	At least one of the two targets <b>Pos</b>		Contamination with a negative or a positive sample during PCR plate preparation? or Master Mix / Water contamination?
<b>EPC</b> <i>D. viviparus</i> PCR external positive control  <b>MANDATORY</b> <i>IN ABSENCE OF MRI</i>	Pos*	Neg	Valid
	Neg	Neg	Problem during qPCR plate preparation: Master Mix error? EPC omission?
	Pos*	Pos	Contamination with a sample during qPCR plate preparation?
<b>Sample process positive Control MRI</b>  <b>RECOMMENDED</b> <i>IF AVAILABLE</i>	Pos <sup>†</sup>	Pos <sup>‡</sup>	Valid
	Neg	Neg	Problem during qPCR plate preparation: Master Mix error? Nucleic acids extract omission or extract not in contact with Master Mix? Process drift: extraction and/or qPCR ?
	Neg	Pos <sup>‡</sup>	Process drift: extraction (in case of exogenous IPC addition directly into qPCR plate and not during extraction) Problem with MRI preparation? Degradation of the sample process positive control?

\* The Ct value obtained must be conform with the value indicated on the Certificate of Analysis (CA).

† The Ct value must be included within control card limits.

‡ The Ct value obtained depends on the thermal cycler and on the extraction protocol used. Ct values for IPC using the recommended extraction protocols are available upon request. BioSella recommends you to determine your own maximal IPC Ct value depending on your own extraction method and thermocycler.

## Samples Reading

**Table 8. Different types of results obtained for the samples**

Targets		Interpretation
<i>D. viviparus</i> (FAM)	Exogenous IPC (Cy5)	
Neg	Pos*	Negative or Undetected
Pos		Positive or Detected
Pos	Neg or Ct>35	<b>Positive or Detected</b> Problem during the IPC addition? Presence of inhibitors †? Competition with the main target?
Neg	Neg or Ct>35	<b>Uninterpretable = Repeat the analysis</b> Nucleic acids extract omission or extract not in contact with Master Mix? Presence of inhibitors †? Nucleic acids degradation in the sample? Problem during the IPC addition? Extraction problem?

\* The obtained Ct value depends on the thermal cycler and the used extraction protocol. This value must be, at least, included within the specified range in the certificate of analysis (CA). Ct values for IPC using the recommended extraction protocols are available upon request. BioSella recommends you determine your own maximal IPC Ct value depending on your own extraction method and thermal cycler.

† In case of inhibition suspicion, 1) Repeat the qPCR with the dilution of extracted nucleic acids at 1/10 or 1/100 in the DNase/RNase free water. 2) Restart the analysis from the extraction step.

Notes :



[www.biosellal.com](http://www.biosellal.com)

### Support Technique

[tech@biosellal.com](mailto:tech@biosellal.com)

+33 (0) 4 26 78 47 62

### Renseignements et commandes

[contact@biosellal.com](mailto:contact@biosellal.com)

+33 (0) 4 26 78 47 60

