

## HANDBOOK

# Bio-T kit® BVDV & BHV1gE

Cat. N° BIOTK106 - 50 reactions

**Detection of Bovine Viral Diarrhea Virus (BVDV) and of Bovine Herpes Virus type 1 (BHV1gE)  
by real-time RT-PCR (qRT-PCR)  
with Endogenous internal positive control (IPC)**

## BOVINE

### Sample types

- Trans-tracheal aspiration liquid (TTA)
- Deep Nasopharyngeal swab (DNS)
- Alveolar bronchial washing (ABW)
- Organs (lungs)
- Individual analysis or by pool up to 3 according to the matrix

### Recommended nucleic acids (NA) extractions

- Magnetic beads extraction (e.g.: BioSella – BioExtract® SuperBall® Cat. N° BES384)
- Silica membrane columns extraction (e.g.: BioSella – BioExtract® Column Cat. N° BEC050 or BEC250)

*Veterinary use only*





## DOCUMENTS MANAGEMENT

The Bio-T kit® BVDV & BHV1gE has two technical handbooks:

- The extraction handbook shared between all the Bio-T kit® of the RESPIRATORY line, displaying BioSella's recommended extraction protocols for each type of sample.
- The Bio-T kit® BVDV & BHV1gE qRT-PCR handbook, presenting the instruction information to perform the qRT-PCR.

The last versions in use for each handbook are indicated on the certificate of analysis (CA) provided with the Bio-T kit® BVDV & BHV1gE.

## MODIFICATIONS MANAGEMENT

BioSella indicates modifications done to this document by highlighting them using the rules presented in the Table below:

MODIFICATIONS MANAGEMENT			
Type of modification Highlighting color	Minor modifications	Type 1 Major modifications	Type 2 Major modifications
Impact on revision / version	Change of revision date No change of version	Change of revision date + change of version	Change of revision date + change of version
Examples of modifications	Corrections: typographical, grammatical or turns of phrase	EPC reference modification	Modification of Master Mix composition
	Addition of new sample type for extraction	Exogenous IPC reference modification	Modification of validated extraction protocol
	Addition of information giving more details or alternative protocol		
	Addition/Suppression of optional information		

## PRESENTATION

### Recommendations for sampling, shipping and storage of samples

Real-time RT-PCR is a powerful technique allowing the detection of few amounts of pathogen genome. Genome can be rapidly degraded depending on the pathogen nature (bacteria / parasites, enveloped viruses...), the genome nature (DNA / RNA) and the sample type (presence of DNase / RNase). Thus, BioSella recommends the following instructions to guarantee an optimal diagnosis.

#### Sampling

To prevent cross-contamination between samples leading to false positive results, it is mandatory to use disposable materials for single use and to avoid direct contact between specimens.

#### Shipping

It is mandatory to ship immediately after sampling or by default to store it at  $\leq -16^{\circ}\text{C}$ . Shipment has to be done within 24h under cover of positive cold.

#### Storage after reception

It is recommended to immediately analyze samples after receipt or freezing at  $\leq -16^{\circ}\text{C}$  for a few months and  $\leq -65^{\circ}\text{C}$  beyond 1 year.

## RESPIRATORY Line

This kit belongs to the RESPIRATORY line which gather a set of kits sharing common extraction and qRT-PCR protocols. It is compatible with BioSella's other kits except with the ones belonging to the PIG and AVIAN lines. (information available on [www.biosellal.com](http://www.biosellal.com)).

In addition to the kits belonging to the RESPIRATORY line, BioSella offers other Bio-T kit® and/or BioLisa® kits allowing the diagnostic of other pathogens involved in Bovine respiratory disease (BRD) such as BVDV or BHV1 individually. For more informations, contact BioSella ([contact@biosellal.com](mailto:contact@biosellal.com))

## Description of the Bio-T kit® BVDV & BHV1gE

The **Bio-T kit® BVDV & BHV1gE** (Cat. N° BIOTK106) contains a ready to use **one-step RT-PCR Master Mix** allowing the detection **in the same reaction well of**:

- **Bovine Viral Diarrhea Virus (BVDV)** with a 6-FAM labelling,
- **Bovine Herpes Virus type 1** non deleted on gE gene (BHV1gE) with a VIC labelling,
- An **Endogenous internal positive control IPC** (gapdh), with a Cy5 labelling, to assess the presence of sufficient amount of host cells, sample integrity, nucleic acids extraction quality and absence of RT- PCR inhibitors.

**Table 1. Description of target viruses**


Cible	BVDV (Fam)	BHV1gE (Vic)*
<b>BHV1 Wild Type Strain</b>	Undetected	Detected
<b>BHV1 vaccine gE deleted</b>	Undetected	Undetected
<b>BHV1 vaccine not deleted on gE</b>	Undetected	Detected
<b>BVDV and other ruminants pestiviruses</b>	Detected	Undetected
<b>Killed BVDV vaccine</b>	Undetected	Undetected
<b>Live BVDV vaccine</b>	Detected	Undetected

\*Cervid Herpes virus type 2 could also be slightly detected by the BHV1gE detection system. However, it is unlikely that cvHV2 causes BRD. In case of doubt, contact the technical support of BioSella ([tech@biosella.com](mailto:tech@biosella.com)).

This kit is based on qualitative detection (detected or non detected) from Trans-tracheal aspiration liquid (TTA), Deep Nasopharyngeal swab (DNS), Alveolar bronchial washing (ABW) and Organs (lungs).

**Extraction protocols recommended by BioSella are described in the extraction handbook of the RESPIRATORY line.**

## Description of the whole process



RESPIRATORY line Extraction handbook		qRT-PCR handbook of the Bio-T kit® BVDV & BHV1gE		
Trans-tracheal aspiration liquid (TTA)*  Deep Nasopharyngeal swab (DNS)*  Alveolar bronchial washing (ABW)*  Organs (lungs)*	BioExtract® SuperBall®  BioExtract® Column	Ready-to-use Master Mix MMBVDBHV1-A	Samples  NC/NCS  Process Positive Control  EPC (EPCBVDBHV1-A)	Dyes: FAM/VIC/Cy5  Passive reference: ROX  Program: Classical program with RT Standard ramping

\* pretreatment mandatory

## Kit contents and storage

**Table 2. Description of the kit contents**

Description	Reference	Volume/tube	Presentation	Storage
<b>Master Mix (MM)</b> Ready to use	MMBVDBHV1-A	750 µl	Transparent cap tube Bag A	≤-16°C Protected from light, « MIX » Area
<b>External Positive Control (EPC)</b> Positive PCR control for BHV1gE and BVDV	EPCBVDBHV1-A	110 µl	Red cap tube Bag B	≤-16°C « Addition of Nucleic acids » Area
<b>Water</b> RNase/DNase free	Aqua-A	1 ml	Blue cap tube Bag B	5°C ± 3 or ≤-16°C « Addition of Nucleic acids » Area

Kit reagents are stable until the expiration date stated on the label, subject to compliance with good storage conditions

## List of consumables and reagents not included in kit

**Table 3. Consumables and reagents not included in kit**

Consumables/ Reagents	Description	Provider	Cat. N°
<b>BioExtract® Column</b>	DNA/RNA column extraction kit (50)	BioSellal	BEC050
<b>BioExtract® Column</b>	DNA/RNA column extraction kit (250)	BioSellal	BEC250
<b>BioExtract® SuperBall®</b>	DNA/RNA Magnetic beads extraction kit (4 x 96)	BioSellal	BES384

For consumables related to the thermal cycler, refer to the user manual of the device.

## Main critical points

- Wear appropriate personal protective equipment (lab coat, disposable gloves frequently changed).
- Work in dedicated and separate areas to avoid contamination: "Extraction" (unextracted samples storage, extraction equipment area), "Mix" (ready to use MM storage, qRT-PCR plates preparation), "Nucleic acids Addition" (Nucleic Acids storage and addition of extracted nucleic acids and controls in the qRT-PCR plate), "PCR" (final area containing the thermal cycler(s)).
- Use dedicated equipment for each working area (gloves, lab coat, pipettes, vortex, ...).
- Use filter tips.
- Before use, thaw all components at room temperature.
- **One-step RT-PCR Master-Mix is less stable than PCR Master-Mix. To guarantee its optimal performance, it is mandatory to extemporaneously defrost the tubes just before the use, to vortex it, to keep it at  $5^{\circ}\text{C} \pm 3$  during the deposit and to refreeze it immediately afterwards.**
- Vortex and spin briefly (mini-centrifuge) all reagents before use.
- Avoid the repetition of freezing-thawing cycles for samples, lysates, extracted nucleic acids.
- **Pathogen's genome detected by the RESPIRATORY line's kits can be DNA or RNA. Working with RNA is more demanding than working with DNA** (RNA instability and omnipresence of the RNases). For these reasons, special precautions must be taken:
  - o Always wear gloves, change them frequently, especially after contact with skin or work surfaces.
  - o Treat all surfaces and equipment with RNases inactivation agents (available commercially).
  - o When wearing gloves and after material decontamination, minimize the contact with surfaces and equipment in order to avoid the reintroduction of RNases.
  - o Use "RNase free" consumable.
  - o It is recommended to store the RNA at  $\leq 5^{\circ}\text{C} \pm 3$  during the manipulation and then freeze it as soon as possible, preferably at  $\leq -65^{\circ}\text{C}$  or by default at  $\leq -16^{\circ}\text{C}$ .
  - o Open and close tubes one by one in order to limit the opening times and avoid any contact with RNases present in the environment (skin, dust, working surfaces...).



# DETECTION OF BVDV & BHV1gE BY qRT-PCR WITH BIOTK106 kit

## Global Procedure

### 1) Establish qRT-PCR plate setup defining each sample position and including the following controls:

- **Negative Control Sample (NCS):** water (or PBS) replaces the sample from the first step of sample preparation.  
This control is mandatory for each extraction series.
- **Negative Amplification Control (NC):** 5 µl of water RNase/DNase free (Aqua-A tube, **blue** cap) replaces sample Nucleic Acids extract on qRT-PCR plate.  
This control is recommended when using the kit for the first time or to verify the absence of Master Mix contamination.
- **External Positive Control of BHV1gE and BVDV (EPC):** Synthetic DNA (tube **EPCBVDV-BHV1-A**, **red** cap), containing specific target for BVDV and BHV1gE.  
This control is mandatory.

**⚠ CAUTION:** *EPC tube handling represents nucleic acids contamination hazard, it is thus recommended to open and handle it in a restricted area, away from other PCR components and to take precautions to avoid cross-contamination with nucleic acids extracts during deposit on the qRT-PCR plate.*

- If available, a **Process Positive Control (MRI)**, a weak positive sample is extracted in parallel with tested samples. After qRT-PCR, MRI Ct values will be monitored on a Shewhart control card. Obtaining conform Ct values validates the whole process. In this case, the use of the EPC, provided with the kit, is not mandatory.

### 2) qRT-PCR plate preparation

In the "MIX" dedicated area

1. After thawing, vortex and rapid centrifugation, **transfer 15 µl Master Mix MMBVDV-BHV1-A** (transparent cap) in each well of interest (samples and controls).

**⚠ NOTE:** *One-step RT-PCR Master-Mix is less stable than PCR Master-Mix. To guarantee its optimal performance, it is mandatory to extemporaneously defrost the tubes just before the use, to vortex it, to keep it at 5°C ± 3 during the deposit and to refreeze it immediately afterwards.*

In the “Nucleic Acids addition” dedicated area

2. **Add 5 µl of extracted nucleic acids (or NCS, water, MRI, EPC: EPCBVDHBV1-A red cap tube)** in each well of interest. Make sure to pipet out in the bottom of the well, in the Master Mix, and to avoid the formation of bubbles.
3. **Seal the plate with an optically clear sealer or close the strip caps.**

In the “PCR” amplification dedicated area

4. **Define the thermal cycler parameters** (see Table 4, Table 5, Table 6 and Table 7)
5. It is recommended to **spin the plate down prior to place it in the thermal cycler**, to prevent drops in the well pit walls.
6. Start the qRT-PCR program. Approximate run time: 80min.

### 3) Thermal cycler settings

This kit was developed and validated on AriaMx™ (Agilent Technologies, Fast ramping by default) and confirmed on ABI PRISM® 7500 Fast (Applied Biosystems) in standard ramping. It is compatible with all thermal cyclers with at least 6-FAM, VIC and Cy5 channels. For more information, contact our technical support.

Table 4. Thermal cycler configuration		
	ABI PRISM® 7500 Fast	AriaMx™
Mode	Quantitation – Standard curve	Quantitative PCR, Fluorescence Probe
Ramping	Standard Ramping	Ramping Fast by default
Passive Reference	ROX	ROX

Table 5. Thermal cycler Settings			
Target	Detectors		Final Volume / well
	Reporter	Quencher	
BVDV	FAM	NFQ-MGB or None*	20 µl  = 15 µl Master Mix + 5 µl extracted nucleic acids or controls <sup>†</sup>
BHV1gE	VIC	NFQ-MGB or None*	
Endogenous IPC	Cy5	NFQ-MGB or None*	
To assign to samples and controls <sup>†</sup>			

\* Depends on the thermal cycler model. Do not hesitate to contact the BioSellal Technical Support (tech@biosellal.com)

† Controls are NC (water), NCS (extracted water) and EPC and or extracted MRI.

**Table 6. CLASSICAL Amplification program settings with RT**

Standard ramping		
Cycles	Time	Temperature
1 cycle	20 min	50°C
1 cycle	5 min	95°C
40 cycles	15 sec	95°C
	30 sec*	60°C
	+ data acquisition	

\* Set 31s for some thermal cyclers such as ABI PRISM® 7500.

NB: This amplification program is compatible with all Bio-T kit\* except for ones belonging to the PIG and AVIAN LINES.

For thermal cycler such as LightCycler®480 and LightCycler®96 (Roche Life Science), it is recommended to use the following program:

**Table 7. PIG/AVIAN Amplification program settings with RT**

Standard ramping		
Cycles	Time	Temperature
1 cycle	20 min	50°C
1 cycle	5 min	95°C
40 cycles	10 sec	95°C
	45 sec	60°C
	+ data acquisition	

NB: This amplification program is compatible with all Bio-T kit\* of the PIG and AVIAN LINES.

## RESULTS INTERPRETATION

To analyze and interpret the signals obtained by qRT-PCR, the Threshold must be set up.

The threshold must be assigned carefully in order to obtain the most reproducible result between different manipulations according to the requirements defined in Annex C of the French Standard **NF U47-600 (part 1)**. A consistent set of positives controls, usually an In-house Reference Material (MRI) or the EPC, is used to set the threshold value above the baseline and in the exponential amplification phase of the plot.

The Threshold Cycle, named « Ct » or « Cq » (depending on thermal cyclers), corresponds to the intersection between the amplification curves and the threshold line. It allows the relative measurement of the concentration of the target in the PCR reaction when a calibrated extract is analyzed in the same series.

The qRT-PCR series is validated if the controls (EPC, MRI, NCS and NC) present valid results, then the result of each sample can be interpreted.

## Main Scenarios

### Controls Reading

Table 8. PCR Controls results interpretation

	Targets			Interpretation
	BVDV (FAM)	BHV1gE (VIC)	Endogenous IPC (Cy5)	
<b>NCS</b> Negative Control Sample <b>MANDATORY</b>	Neg	Neg	Neg	Valid
	At least one of the three targets <b>Pos</b>			Contamination with a positive/negative sample during extraction step or during qRT-PCR plate preparation.
<b>NC</b> Negative PCR Control <b>OPTIONAL</b>	Neg	Neg	Neg	Valid
	At least one of the three targets <b>Pos</b>			Contamination with a positive/negative sample during extraction step or during qRT-PCR plate preparation or Master Mix/water contamination.
<b>EPC</b> BVDV and BHV1gE PCR external positive control <b>MANDATORY</b>	Pos*	Pos*	Neg	Valid
	Neg	Neg	Neg	Problem during qRT-PCR plate preparation: Master Mix error? EPC omission?
	Pos*	Pos*	Pos	Contamination with a sample during qRT-PCR plate preparation?
<b>Sample process positive Control MRI</b> <b>RECOMMENDED</b> <i>IF AVAILABLE</i>	Pos <sup>†</sup>	Pos <sup>†</sup>	Pos <sup>‡</sup>	Valid
	Neg	Neg	Neg	Problem during qRT-PCR plate preparation: Master Mix error? Nucleic acids extract omission or extract not in contact with Master Mix? Process drift: extraction and/or qRT-PCR ? Degradation of the sample process positive control?
	Neg	Neg	Pos <sup>‡</sup>	Process drift: extraction Problem with MRI preparation? Degradation of the sample process positive control?

\* The Ct value obtained must be conform with the value indicated on the Certificate of Analysis (CA).

<sup>†</sup> The Ct value must be included within control card limits.

<sup>‡</sup> The obtained Ct value depends on the thermal cycler, the sample type and the used extraction protocol. IPC Ct values for recommended extraction protocols are available upon request. BioSella recommends you determine your own maximal IPC Ct value depending on your own extraction method and thermal cycler.

Note:

Endogenous IPC targets a gene expressed by bovine cells, thus it cannot be detected in NCS, NC and EPC. However, due to cross-reaction between ruminant GAPDH and human GAPDH, a slight signal can be observed for IPC in the controls, the Ct value of this signal must be lower than 35.

## Samples Reading

**Table 9. Different types of results obtained for the samples**

Targets			
BVDV (FAM)	BHV1gE (VIC)	Endogenous IPC (Cy5)	Interpretation
Neg	Neg	Pos*	Negative or Undetected
Pos	Pos		Positive or Detected
At least one of two targets <b>Pos</b>			Positive or Detected for the positive target Negative or Undetected for the negative target
Pos	Pos	Neg or Ct>35	Positive or Detected Lack of host cells? Presence of inhibitors¹? Competition with the main target?
At least one of the target is <b>Neg</b>		Neg or Ct>35	Positive or Detected for the positive target Uninterpretable for the negative target= Repeat the analysis. Presence of inhibitors¹? Nucleic acids degradation in the sample? Sampling problem: lack of cells? Extraction problem?
Neg	Neg	Neg or Ct>35	Uninterpretable = Repeat the analysis Problem during qRT-PCR plate preparation: Master Mix error? Nucleic acids extract omission or extract not in contact with Master Mix? Presence of inhibitors¹? Nucleic acids degradation in the sample? Sampling problem: lack of cells? Extraction problem?

\*The obtained Ct value depends on the thermal cycler, the sample type and the used extraction protocol. IPC Ct values for recommended extraction protocols are available upon request. BioSella recommends you determine your own maximal IPC Ct value depending on your own extraction method and thermal cycler.

† In case of inhibition suspicion, 1) Repeat the qRT-PCR with the dilution of extracted nucleic acids at 1/10 or 1/100 in the DNase/RNase free water. 2) Restart the analysis from the extraction step.



**[www.biosellal.com](http://www.biosellal.com)**

### **Technical Support**

**[tech@biosellal.com](mailto:tech@biosellal.com)**

**+33 (0) 4 26 78 47 62**

### **Information and orders**

**[contact@biosellal.com](mailto:contact@biosellal.com)**

**+33 (0) 4 26 78 47 60**

