SVANOVIR® TGEV/PRCV-Ab

**Transmissible Gastroenteritis Virus/
Porcine Respiratory Coronavirus
Differentiating Test**
Antibody Test

**Transmissibles Gastroenteritis Virus/
Porcines Respiratorisches Coronavirus Differenzierungstest**

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Contents	Art.No. 10-7500-02
Microtitre plate Microtitre plates (96 wells) coated with non-infectious TGEV antigen (sealed and stored dry)	2 (Strips) 12 x 8
Conjugate Concentrate (horseradish peroxidase conjugated anti-mouse IgG antibodies)	1 x 30 µL
PBS-Tween Solution 20x concentrate	1 x 125 mL
Conjugate Dilution Buffer - Contains preservatives	1 x 25 mL
Substrate Solution (Tetramethylbenzidine in substrate buffer containing H ₂ O ₂) - STORE IN THE DARK!	1 x 20 mL
Stop Solution Contains sulphuric acid (2M) - CORROSIVE!	1 x 10 mL
Anti-TGEV mAb Solution - Contains 0.07% ProClin	1 x 10 mL
Anti-TGEV/PRCV mAb Solution - Contains 0.07% ProClin	1 x 10 mL
A. Positive TGEV Control Serum - Contains 0.07% ProClin	1 x 2.5 mL
B. Negative TGEV/PRCV Control Serum - Contains 0.07% ProClin	1 x 2.5 mL
C. Positive PRCV Control Serum - Contains 0.07% ProClin	1 x 2.5 mL

This manual covers the following
SVANOVIR® TGEV/PRCV-Ab kit:
Article number 10-7500-02

Transmissible Gastroenteritis Virus/Porcine Respiratory Coronavirus Differentiating Test Antibody Test

Name and Application

The SVANOVIR® TGEV/PRCV-Ab test kit is an Enzyme Linked Immunosorbent Assay (ELISA) designed to differentiate between TGEV and PRCV specific antibodies in porcine serum. A maximum of 88 porcine sera can be tested with the kit.

General information

Transmissible gastroenteritis (TGE) is a highly contagious viral enteric disease of swine characterised by vomiting, watery diarrhoea, anorexia and dehydration. Mortality is high in younger piglets. Suckling pigs older than 3 weeks usually survive the disease. The clinical signs in adult pigs are limited to loss of appetite and diarrhoea. The disease is caused by the coronavirus TGEV (transmissible gastro-enteritis virus). Ingested TGEV replicates primarily in the villous epithelial cells of the small intestine¹. Milk and faeces are the main sources of virus transmission from carrier animals. Starlings and houseflies can act as mechanical vectors for the virus. The porcine respiratory coronavirus (PRCV) is antigenically related to TGEV but replicates primarily in the tonsils and the respiratory tract^{2,3}. PRCV may cause pneumonia in swine⁴. Normally PRCV induce only mild or no clinical respiratory signs^{5,6}. Because PRCV and TGEV have many similar antigenic sites on the surface of the virus^{7,8} and the problems encountered with carrier animals, it is important to clearly differentiate between these two viruses in a routine diagnostic test.

Principle

The kit procedure is based on the blocking ELISA. In this procedure samples are exposed to non-infectious TGEV antigen-coated wells in microtitre plates or strips. If anti-TGEV or anti-PRCV antibodies are present in the test sample they will bind to viral antigens in the well and block these antigenic sites. If anti-TGEV or anti-PRCV antibodies are absent in the sample these sites will remain free. When anti-TGEV or anti-TGEV/PRCV monoclonal antibodies are added they will attach to specific free sites on the virus. By adding horseradish peroxidase (HRP) conjugated anti-mouse Ig antibodies they will bind to the monoclonal antibodies present in the well. When the substrate is added a blue colour develops. A negative result is indicated by a strong colour change. The reaction is stopped by addition of stop solution; the colour changes to yellow. The optical density (OD) is measured at 450 nm in a microplate photometer and the percent inhibition of the positive controls (A, C) and test samples relative to the negative control (B) is calculated.

Materials needed but not provided

1. Precision pipettes
2. Disposable pipette tips
3. Distilled, deionised or any similar high quality water
4. Wash bottle, multichannel pipettor or plate washer
5. Containers: 1 to 2 litres for PBS-Tween and 20 mL for Conjugate Solution
6. Microplate photometer, 450 nm filter

Specimen information

2 x 50 µL of undiluted blood serum is required for each sample. Each test sample requires the use of two sample wells. Fresh, refrigerated or previously frozen serum can be used.

Preparation of reagents

PBS-Tween Buffer:

Dilute the PBS-Tween Solution 20 x concentrate 1/20 in distilled water. Prepare 500 mL per plate by adding 25 mL PBS-Tween solution to 475 mL distilled water and mix thoroughly.

N.B. Please check that there is no crystal precipitation in the bottle. If crystals are seen, please warm and shake well.

HRP Conjugate Solution:

Dilute the HRP Conjugate Concentrate in Conjugate Dilution Buffer according to the given dilution factor (f) on the label of the conjugate vial. Prepare the reagent immediately before use and mix it thoroughly on a shaker. Prepare 10 mL for one 96 well plate or 1 mL for each 8 well strip. Discard any unused solution.

Precautions

1. Carefully read and follow all instructions.
2. Store the kit and all reagents at 2-8°C (36-46°F).
3. All reagents should equilibrate to room temperature 18-25°C (64-77°F) before use.
4. Handle all materials according to the Good Laboratory Practice.
5. Do not mix components or instruction manuals from different test kit batches.
6. Care should be taken to prevent contamination of kit components.
7. Do not use test kit beyond date of expiry.
8. Do not eat, drink, or smoke where specimens or kit reagents are handled.
9. Use a separate pipette tip for each sample.
10. Do not pipette by mouth.
11. Include positive and negative control sera (3 controls) on each plate or test strip series.
12. Use only distilled deionised or any similar high quality water for preparation of reagents.
13. When preparing the buffers, etc., measure the required volume.
14. The Stop Solution contains sulphuric acid, which is corrosive.
15. All unused biological materials should be disposed according to the local, regional and national regulations.

Recommendations!

The volume of reagents is sufficient for at least 10 separate test occasions. Strips with broken seal can be stored at 2-8°C (36-46°F) for up to 4 weeks.

Procedure

1. All reagents should equilibrate to room temperature 18-25°C (64-77°F) before use. Label each strip with a number.

2. Wash the plates/strips 3 times with PBS-Tween buffer: fill up the wells at each wash, leave the plate with the buffer 1 minute, empty the plates and tap hard to remove all remains of fluid.

3. Add 100 µL of Positive Control Serum (A and C) and 100 µL of Negative Control Serum (B) to selected wells. It is recommended to run the control sera in duplicates — see figure 1.

4. Add 50 µL of PBS-Tween and 50 µL of each undiluted serum sample to selected wells. Add samples in duplicates.

5. Shake the plate thoroughly by tapping it on the side. Seal the plate/strips and incubate for 2 hours at +37°C (98.6 °F).

6. Repeat step #2.

7. Add 100 µL Anti-TGEV Solution to one of the test sample wells and add 100 µL Anti-TGEV/PRCV Solution to the other test sample well. Repeat this step for each pair of controls and samples (see figure #1).
8. Seal the plate/strips and incubate for 30 minutes at room temperature 18-25°C (64-77°F).

9. Repeat step #2.

10. Add 100 µL of diluted HRP Conjugate Solution to each well, seal the plate/ strips and incubate for 30 minutes at room temperature 18-25°C (64-77°F).

11. Repeat step #2.

12. Add 100 µL Substrate Solution to each well and incubate for 30 minutes at room temperature 18-25°C (64-77°F). Begin timing when the first well is filled.

13. Stop the reaction by adding 50 µL Stop Solution to each well and mix content of wells thoroughly. Add the Stop Solution in the same order as the Substrate Solution in step #12.

14. Measure the optical density (OD) of the controls and samples at 450 nm in a microplate photometer (use air as blank). Measure the OD within 15 minutes after the addition of Stop Solution to prevent fluctuations in OD values.

Figure #1: Recommended process chart for addition of controls and samples

	1	2	3	4
A	Control A	Control A	Sample 5	Sample 5
B	Control B	Control B	Sample 6	Sample 6
C	Control B	Control B	Sample 7	Sample 7
D	Control C	Control C	Sample 8	Sample 8
E	Sample 1	Sample 1	Sample 9	Sample 9
F	Sample 2	Sample 2	Sample 10	Sample 10
G	Sample 3	Sample 3	etc.	etc.
H	Sample 4	Sample 4	etc.	etc.
	Anti-TGEV solution	Anti-TGEV/PRCV solution	Anti-TGEV solution	Anti-TGEV/PRCV solution

Calculations

N.B. All calculations must be done for the TGEV and the TGEV/PRCV mAb solution separately.

1. Calculate the mean optical density (OD) value for each of the controls and samples.
2. Calculate the percent inhibition (PI) values for the positive controls as well as samples, using the following formula:

$$PI = \frac{OD_{Neg. control} - OD_{Sample or Pos. control}}{OD_{Negative control}} \times 100$$

Interpretation of the results

Criteria for test validity

To ensure test validity the TGEV/PRCV Negative Serum Control (B) should have an OD value greater than 0.5 and the Control Sera A and C as follows:

% Inhibition (PI)		
	Anti-TGEV	Anti TGEV/PRCV
A	> 60	> 60
C	< 15	> 60

Interpretation

PI	> 60	positive	(+)
PI	45-60	doubtful	(+/-)
PI	< 45	negative	(-)

Serum samples with any PI-values falling in the doubtful zone(+/-) should be interpreted as inconclusive or invalid and assay must be repeated. If still inconclusive or invalid a new sample with at least 3 weeks interval should be tested.

Anti-TGEV	Anti-TGEV/PRCV	Interpretation
-	-	Negative for both TGEV and PRCV
-	+	Positive for PRCV (negative for TGEV)
+	+	Positive for TGEV
+/-	-	Invalid test result*
+	-	Invalid test result*
+	+/-	Invalid result**
+/-	+	Inconclusive for TGEV
-	+/-	Inconclusiv for PRCV
+/-	+/-	Inconclusive for TGEV and PRCV

* Some sera with high levels of antibodies to PRCV may show a negative reaction in the blocking ELISA due to a prozone phenomena. These sera will score positive when tested at a higher dilution (for example 1/250 or 1/500).

**Invalid test result due to either prozone phenomena (see*), or technical error in performing the test.

Data from Canadian study (Herd basis)

Study Reference	Sample Size (nos. of herds) ¹ Kit/reference ²				Sample Type	Relative Sensitivity and Specificity (95% Confidence Interval)
	+/+	-/+	+/-	-/-		
J. Vet. Diag. Invest. (2002) 14:97-105 ⁹	28	2	4	66	Serum	Sensitivity: 93.3 (77.9%, 99.2%) Specificity: 94.3% (86%, 98.4%)

1. Known serological status

2. Clinical signs

References

1. Saif, L. J. and Bohl, E. H. (1986) Transmissible Gastroenteritis. In Diseases of Swine 6th Edition. Edited by A.D. Leman, R.D. Glock, W. C. Mengeling, R.H. Penny, E. Scholl, and B. Straw. Iowa State Uni. Press, Ames, Iowa U.S.A. pp. 255-274.
2. Furuuchi, S., Shimizu, Y., and Kumagai, T. (1978/1979) Multiplication of low and high cell culture passaged strains of transmissible gastroenteritis virus in organs of newborn piglets. Vet. Microb. 3, 169-178.
3. Kemeny, L. J., Wiltsey, V. L., and Riley, J. L. (1975) Upper respiratory infection of lactating sows with transmissible gastroenteritis virus following contact exposure to infected pig-lets. Cornell Vet. 65, 352-362.
4. van Nieuwstadt, A. P. and Pol, J. M. A. (1989) Isolation of a TGE virus-related respiratory coronavirus causing fatal pneumonia in pigs. Vet. Rec. 124, 43-44.
5. Brown, I. and Cartwright, S. F. (1986) New porcine coronavirus? Vet. Rec. 119, p. 282.
6. Pensaert, M., Callebaut, P. and Vergote, J. (1986) Isolation of a porcine respiratory non-enteric coronavirus related to transmissible gastroenteritis. The Veterinary Quarterly. 8, 257-260.
7. Garwes, D. J., Stewart, F., Cartwright, S. F. and Brown, I. (1988). Differentiation of porcine coronavirus from transmissible gastroenteritis virus. Vet. Rec. 122, 86-87.
8. Callebaut, P. *et al.* (1988). Antigenic differentiation between TGEV of swine and a related porcine respiratory coronavirus. J. Gen. Virol.
9. Carman, S. *et al.* (2002). Field validation of a commercial blocking ELISA to differentiate antibody to transmissible gastroenteritis virus (TGEV) and porcine respiratory coronavirus and to identify TGEV-infected swine herds. J. Vet. Diagn. Inves.t 14:97-105.

Produced by SVANOVA
Distributed by BIOSELLAL SAS
Bâtiment Accinov
317 Avenue Jean Jaurès
69007 LYON
FRANCE
+33 (0) 4 26 78 47 60

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Referenzen

1. Saif, L. J. and Bohl, E. H. (1986) Transmissible Gastroenteritis. In Diseases of Swine 6th Edition. Edited by A.D. Leman, R.D. Glock, W. C. Mengeling, R.H. Penny, E. Scholl, and B. Straw. Iowa State Uni. Press, Ames, Iowa U.S.A. pp. 255-274.
2. Furuuuchi, S., Shimizu, Y., and Kumagai, T. (1978/1979) Multiplication of low and high cell culture passaged strains of transmissible gastroenteritis virus in organs of newborn piglets. Vet. Microb. 3, 169-178.
3. Kemeny, L. J., Wiltsey, V. L., and Riley, J. L. (1975) Upper respiratory infection of lactating sows with transmissible gastroenteritis virus following contact exposure to infected pig-lets. Cornell Vet. 65, 352-362.
4. van Nieuwstadt, A. P. and Pol, J. M. A. (1989) Isolation of a TGE virus-related respiratory coronavirus causing fatal pneumonia in pigs. Vet. Rec. 124, 43-44.
5. Brown, I. and Cartwright, S. F. (1986) New porcine coronavirus? Vet. Rec. 119, p. 282.
6. Pensaert, M., Callebaut, P. and Vergote, J. (1986) Isolation of a porcine respiratory non-enteric coronavirus related to transmissible gastroenteritis. The Veterinary Quarterly. 8, 257-260.
7. Garwes, D. J., Stewart, F., Cartwright, S. F. and Brown, I. (1988). Differentiation of porcine coronavirus from transmissible gastroenteritis virus. Vet. Rec. 122, 86-87.
8. Callebaut, P. *et al.* (1988). Antigenic differentiation between TGEV of swine and a related porcine respiratory coronavirus. J. Gen. Virol.
9. Carman, S. *et al.* (2002). Field validation of a commercial blocking ELISA to differentiate antibody to transmissible gastroenteritis virus (TGEV) and porcine respiratory coronavirus and to identify TGEV-infected swine herds. J. Vet. Diagn. Inves.t 14:97-105.

Hergestellt von SVANOVA
Verteil von BIOSELLAL SAS
Bâtiment Accinov
317 Avenue Jean Jaurès
69007 LYON

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2. Signos clínicos

Referencias

1. Saif, L. J. and Bohl, E. H. (1986) Transmissible Gastroenteritis. In Diseases of Swine 6th Edition. Edited by A.D. Leman, R.D. Glock, W. C. Mengeling, R.H. Penny, E. Scholl, and B. Straw. Iowa State Uni. Press, Ames, Iowa U.S.A. pp. 255-274.
2. Furuuuchi, S., Shimizu, Y., and Kumagai, T. (1978/1979) Multiplication of low and high cell culture passaged strains of transmissible gastroenteritis virus in organs of newborn piglets. Vet. Microb. 3, 169-178.
3. Kemeny, L. J., Wiltsey, V. L., and Riley, J. L. (1975) Upper respiratory infection of lactating sows with transmissible gastroenteritis virus following contact exposure to infected pig-lets. Cornell Vet. 65, 352-362.
4. van Nieuwstadt, A. P. and Pol, J. M. A. (1989) Isolation of a TGE virus-related respiratory coronavirus causing fatal pneumonia in pigs. Vet. Rec. 124, 43-44.
5. Brown, I. and Cartwright, S. F. (1986) New porcine coronavirus? Vet. Rec. 119, p. 282.
6. Pensaert, M., Callebaut, P. and Vergote, J. (1986) Isolation of a porcine respiratory non-enteric coronavirus related to transmissible gastroenteritis. The Veterinary Quarterly. 8, 257-260.
7. Garwes, D. J., Stewart, F., Cartwright, S. F. and Brown, I. (1988). Differentiation of porcine coronavirus from transmissible gastroenteritis virus. Vet. Rec. 122, 86-87.
8. Callebaut, P. *et al.* (1988). Antigenic differentiation between TGEV of swine and a related porcine respiratory coronavirus. J. Gen. Virol.
9. Carman, S. *et al.* (2002). Field validation of a commercial blocking ELISA to differentiate antibody to transmissible gastroenteritis virus (TGEV) and porcine respiratory coronavirus and to identify TGEV-infected swine herds. J. Vet. Diagn. Inves. 14:97-105.

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Kit distribué par BIOSELLAL SAS