

INFORMATION

SVANOVIR® NEOSPORA-Ab

Cat.No. : 104898

Batches : B73558, B86522, B94821

Change : shorter stability for strips with broken seal

Dear customers :

Following a communication from the manufacturer, we inform you that for the batches listed below, the open plates can be stored at 4°C for 1 week only and not 4 weeks as indicated in the handbook MU/NEO/0002/EN provided in the kits.

For any information, you can contact Biosellal



27 Chemin des Peupliers
69570 Dardilly
FRANCE
+33 (0) 4 26 78 47 60
www.biosellal.com
contact@biosellal.com

